## Instructions For Use

### Nor-NOHA Blood Collection Tube/ Nor-NOHA Blutentnahmeröhrchen/ Nor-NOHA Tubes de Prélèvement Sanguine/ Nor-NOHA de Tubos Extracción de Sangre

In Vitro Diagnostic Use/ Zur in-vitro-diagnostischen Verwendung/ Pour diagnostic in vitro / Para uso diagnóstico in vitro  
REF: 700006

EN

#### Intended Use & Product Description

Nor-NOHA Blood Collection Tube is a non-sterile, non-automated, single use blood collection tube intended for use by healthcare professionals. The Nor-NOHA tube contains an additive mixture including an arginase inhibitor, nor-NOHA, which blocks enzymatic activity of arginase, preventing arginine breakdown during the processing of the whole blood sample to plasma. The Nor-NOHA tube is to be used for arginine analyses in patients treated with medicinal products which exhibit its effect by degrading arginine, to enable accurate analysis of plasma arginine. The tube is not sterile (see Precautions/ Cautions below).

The interior of the tube wall is coated with 18 µg nor-NOHA and 3.6 mg K2EDTA. Nor-NOHA inhibits degradation of arginine and EDTA is an anticoagulant that block the coagulation cascade by binding calcium ions.

Nor-NOHA tubes have a dimension of 13x75 mm, are made of inert plastic (PET), are evacuated and non-sterile. They are for single use and have a draw volume of 2 ml (total capacity 5 ml). The polypropylene cap (white) includes a bromobutyl rubber stopper (grey).

#### Precautions/Cautions

1. Tubes are not sterile. Hence, the tubes should not be used for blood collection via direct phlebotomy unless using a catheter/ butterfly needle of at least 12.5 cm in length, which is equipped with a multi-sample luer adapter (MSLA). This will reduce the possibility of backflush from the non-sterile tube to the sample donor.
2. Do not use tubes if foreign matter is present.
3. If blood is collected through an intravenous (IV) line, ensure that the line has been cleared of IV solution before beginning to fill blood collection tubes. This is critical to avoid erroneous laboratory data from IV fluid contamination.
4. Handle all biological samples according to the policies and procedures of your facility.
5. Biological samples may transmit blood-borne pathogens. Obtain appropriate medical attention in the case of any exposure to biological samples and the resulting infection risk.
6. Discard all accessories in biohazard containers approved for their disposal.
7. The tubes are not intended for sterilisation or autoclaving.
8. Do not re-use.
9. Do not use tubes after their expiration date.

Only applicable for member states of the European Union/ Great Britain: Should any serious incidents occur in relation to the product, these must be reported to the manufacturer and the relevant authority, in which the user/patient is established.

#### Storage, Handling & Disposal

Store tubes at 15-25°C prior to use. Exceeding the maximum recommended storage temperature may lead to impairment of the tube quality. Store in a dry place, away from sources of heat and direct sunlight.

The non-interference of utilizing the tubes on the instrument can be tested and validated directly by the user. Refer to specific documents (i.e., guidelines, standards) for additional information. Follow the analyser manual for the correct handling and ensure methods are validated. The tubes withstand a centrifugation of no more than 3000g for 10 minutes. Disposal: Follow general hygiene guidelines and legal regulations for the proper disposal of infectious material.

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Immedica Pharma AB